

Intelligent Evaporation



MADE IN
GERMANY

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.

Chemical-resistant **vacuum seal** allows **for years of operation** and reduces spare parts budget up to 75 %

Bath power cable coupling complies with protection class IP 67 to **prevent corrosion and short circuits**

Detachable panel allows for the **highest operational safety** from outside closed fume hoods

Continuous unattended evaporation without limits - the automatic module Distimatic Bench-top with automatic drainage of residue is an affordable and more flexible alternative to large-scale rotary evaporators

The flange made of **chemical-resistant PPS** reduces corrosion, maintenance and spare parts cost to a minimum

Non-sticking vapor tube with patented clamping sleeve eliminates broken glass and reduces spare parts cost significantly

The **patented Easy-Clip** provides **an easy removal** of sticking flasks and **can't get lost** due to the direct connection to the drive

Safety bath will power off if temperature overshoots by 5 °C or runs dry

210 °C heating bath with safety handles and pour spout accommodates **5-liter flasks**

A metal support between the heating bath and the base unit **prevents bath instability**

A separate on/off switch for heating **prevents unintentional heat-up** - the button is illuminated for visual control

